

PROJECT CASE STUDY

Custom TruCytes™ Controls Validate Prognostic FcγRIIIa Flow Cytometry Test on Fixed Platelets

CLIENT

Prolocor

THERAPEUTIC AREA

Cardiovascular disease

PROJECT OVERVIEW

Cardiovascular disease remains the leading cause of death and disability worldwide. Antithrombotic treatments such as dual antiplatelet therapy (DAPT) reduces ischemic cardiovascular events by inhibiting platelets, but its use can increase bleeding risk. Physicians must carefully balance ischemia and bleeding risks to make treatment decisions.

Prolocor has a goal to help clinicians assess thrombosis risk and personalize antiplatelet therapy with a novel FcγRIIIa test. FcγRIIIa amplifies platelet activation, and high platelet FcγRIIIa (pFCG) expression correlates with an increased risk of recurrent cardiovascular events, making it a valuable biomarker for patient stratification.

High-risk patients with elevated pFCG levels may benefit from more intensive or prolonged treatment. However, in low-risk patients, the bleeding risks of extended therapy may outweigh its benefits, making a shorter duration the safer, more effective choice.

The Prolocor pFCG™ test quantifies FcγRIIIa levels using a novel 5G1 antibody to bind FcγRIIIa on previously fixed platelets, addressing issues with platelet degradation during in vitro analysis and improving assay precision.¹

However, developing a reliable diagnostic test requires robust baseline assessment and consistent controls. Biological controls present challenges, including donor variability, limited availability, and stability issues, which complicate assay development, validation, and routine clinical implementation.



SLINGSHOT BIOSCIENCES

To overcome these challenges, Prolocor integrated Slingshot Biosciences' TruCytes™ FcγRIIA-conjugated cell mimics as standardized assay controls for their pFCG™ test validation studies. TruCytes™ provided a stable, reproducible alternative to biological specimens (fixed platelets) to improve assay precision, accuracy, and to assess the test's linear range. Implementation of TruCytes™ as internal assay controls improved test consistency across different sample batches and strengthened the reliability of the pFCG™ test across diverse clinical applications.

CUSTOMER CHALLENGES

Achieving assay precision is critical in bioanalysis and diagnostics to generate reliable prognostic results for patient treatment. However, biological variability is a diagnostic challenge—platelets from different donors may express different levels of a given protein or respond differently to stimuli. When researchers rely on biological samples as controls, this variability makes it difficult to establish clear and reproducible benchmarks during assay development. Furthermore, unlike blood, plasma, or serum, platelets cannot be 'spiked' with a given amount of a specific protein. This makes assessment of precision and accuracy more difficult to assess.

SOLUTION

Slingshot Biosciences leverages biochemistry, high-precision manufacturing, and polymer chemistry to engineer cell mimics with tunable size, morphology, and biochemical properties that can be matched to any cell type.

Using this technology, Slingshot produced custom TruCytes™ cell mimics conjugated with purified FcγRIIa at different protein loading levels that span the biologic range of ~500 to 8,000 molecules of FcγRIIa per cell mimic.

This represents a broad range of protein expression values across biological platelet samples from low values seen in healthy subjects to higher values seen in patients with cardiovascular disease and cancer.

An accelerated stability study conducted by Slingshot demonstrated that product stability is maintained for 18 months when stored in lyophilized form at -20°C. Upon reconstitution in phosphate-buffered saline with 0.2% bovine serum albumin, TruCytes™ remained stable for 30 days at 4°C.



RESULTS

Assessment of Assay Precision

pFCG™ test validation with biological specimens (fixed platelets) showed an intra-assay coefficient of variation (CV) of $2.1 \pm 0.1\%$ and inter-assay CVs of $4.5 \pm 1\%$ (intraday) and $6.5 \pm 0.4\%$ (up to 5 days post-fixation) (Figure 1A). In comparison, the Inter-assay CVs were $2\% \pm 0.6\%$ (intraday) and $9.9\% \pm 2.1\%$ (interday) using the TruCytes™ (Figure 1B).




FIGURE 1: Evaluation of pFCG® assay precision. A) The interassay CV assessed across 50 biological samples. Intraday samples were performed on the same day with the second pFCG™ test performed at least 2h after the first test. Interday samples were performed on days 1–5 after fixation of samples. B) The interassay CV assessed with Fc γ RIIIa-conjugated TruCytes™. Intraday samples were performed on the same day with the second pFCG™ test performed at least 2h after the first test. Interday samples were performed on 20 non-consecutive days.

Assessment of Assay Accuracy

Figure 2 demonstrates the use of TruCytes™ to evaluate the accuracy and linearity of Fc γ RIIIa quantification. The Prolocor pFCG™ test demonstrated excellent accuracy and linearity across biologic expression levels with a linear regression of $R^2 = 0.984$.




FIGURE 2: Evaluation of Fc γ RIIIa quantitation and pFCG® assay linearity. A) The pFCG test accurately quantified Fc γ RIIIa protein expression on the TruCytes™. Results were compared with the expected results. B) Linearity of the pFCG™ test was assessed with the TruCytes™. The pFCG™ test was performed (n = 765 determinations) on samples that ranged from 591 to 7640 molecules of Fc γ RIIIa/cell mimic.



SLINGSHOT
BIOSCIENCES

CUSTOMER OUTCOME

The Prolocor pFCG™ test offers rapid and accurate FcγRIIIa quantification on fixed platelets via flow cytometry. Analytical validation studies demonstrate the test's high accuracy, with intra-assay variation at just 2% and inter-assay variation at 10%, both well below the FDA's acceptable threshold of 25%.

The inclusion of FcγRIIIa-conjugated TruCytes™ as internal controls provides a significant advantage over traditional biological controls. TruCytes™ offer a stable and reproducible alternative that outperforms biological specimens, enabling more consistent and precise quantification of FcγRIIIa. This enhanced consistency and reliability make TruCytes™ a key component for maintaining high precision and accuracy in the Prolocor pFCG™ test, further solidifying its value as a prognostic tool in clinical settings.

As clinicians seek better tools for informed decision-making, the prognostic capabilities of the Prolocor pFCG™ test are a critical asset for enabling personalized cardiovascular disease treatment that prioritizes patient safety.

"We want to thank our partners for developing the novel TruCytes™ cell mimics that are central to the test and pivotal to demonstrating the accuracy and power of this test," said Dr. David Schneider, Prolocor co-founder and Chief Science Officer.

REFERENCES

1. Schneider, David J et al. "Assessing prognosis by quantifying FcγRIIIa on fixed platelets." *Bioanalysis* vol. 16,19-20 (2024): 1025-1032. doi:10.1080/17576180.2024.2395706

To learn more about how TruCytes™ can support your diagnostic and bioanalysis needs, visit:

slingshotbio.com/trucytes

