



1. Technical Data Sheet

Summary	SpectraComp® eGFP compensation and calibration controls are state-of-the-art cell mimics with actual eGFP presented on the surface.
Application	<p>SpectraComp® eGFP is intended as a compensation and calibration control for cells expressing eGFP protein. SpectraComp® eGFP kit consists of four distinct populations including a negative and three positive levels. Each positive level varies in the number of eGFP proteins, giving each population a distinct fluorescence intensity level.</p> <p>Note: SpectraComp® eGFP performance has been verified and validated on analytical flow cytometers and not on cell sorters.</p> <p>For Research Use Only. Not for use in diagnostic or therapeutic procedures.</p>
Materials	SpectraComp® eGFP are cell mimics that are suspended in aqueous solution and are packaged in a convenient dropper bottle. Each drop contains approximately 1×10^5 beads.
Handling and Safety	No special handling or safety precautions are necessary. See the Safety Data Sheet (SDS) at www.slingshotbio.com .
Storage	SpectraComp® eGFP should be stored at 2 - 8 °C once the product is received.
Expiration	One year from the date of manufacturing
Instructions for Use	<ol style="list-style-type: none"> 1. Unpack and vortex all vials on high for 2 - 3 seconds to resuspend cell mimics. 2. Label the required number of tubes and add 150 µl of 1X PBS buffer to each tube. 3. Add 1 drop of the negative, low, medium or high beads in the appropriately labeled tubes. See the illustrations below as examples of how the different cell mimic populations can be configured. <p>Do not store cell mimic populations after mixing. Samples are intended to be analyzed immediately after mixing.</p> <p>Configuration 1. One tube containing only negative and one tube</p>








containing only positive(s).



Configuration 2. One tube containing only negative, and one tube containing negative and positive(s).



Configuration 3. One tube containing negative and low positive, one tube containing negative and medium positive, and one tube

<p>containing negative and high positive.</p> 				
<p>4. Vortex to suspend the cell mimics.</p> <p>Note: Protect the samples from light and analyze the samples as soon as possible.</p> <p>5. Acquire the samples using the same Forward and Side Scatter parameters (FSC-A and SSC-A) as would be used for the actual cells. Refer to QC data for expected populations (prepared with configuration 2).</p>				
QC Data	Figure 1. SpectraComp® eGFP Figure 1 (A, B)			
	<table><tr><td></td><td></td></tr><tr><td>Figure 1. (A) Gate on scatter population of SpectraComp® eGFP</td><td>Figure 1. (B) Histogram plot showing negative, low, medium and high signal peaks</td></tr></table>			Figure 1. (A) Gate on scatter population of SpectraComp® eGFP
				
Figure 1. (A) Gate on scatter population of SpectraComp® eGFP	Figure 1. (B) Histogram plot showing negative, low, medium and high signal peaks			
Technical Support	For technical support regarding this product please contact: support@slingshotbio.com			