

# Enhancing Flow Cytometry Assay Performance with Cell Mimics Controls for Immune Cell Profiling and Quality Control



SLINGSHOT  
BIOSCIENCES


Swetha Gunturu, Brian Kerfs, Kanwal Palla, Selina Yee  
Slingshot Biosciences, Inc



## ABSTRACT

Flow cytometry is essential for immune profiling, but a lack of standardized controls limits reproducibility and assay reliability. Cell mimics expressing single, defined biomarkers offer a consistent, flexible alternative to biological controls. These mimics integrate into diverse flow cytometry assays, including Lymphocyte Subset, T-cell activation, exhaustion, CAR-T detection, and post-infusion biomarker analysis. Their stability, reproducibility, and compatibility with multicolor panels enable improved assay calibration, gating validation, and immune monitoring. To evaluate performance, mimics bearing single activation markers (CD25, CD69, CD278, CD38) and single exhaustion markers (PD-1, CTLA-4, TIM-3) were spiked into PBMCs and subjected to standard staining workflows. Resulting percent-positive signals closely aligned with those from activated PBMCs, confirming the mimics' utility as reliable, reproducible assay controls. This platform enables more accurate, standardized T-cell analysis for applications in immunotherapy, vaccine development, and immune monitoring.


## Single Biomarker TruCytes Cell Mimic Controls Can Be Spiked Into Peripheral Blood Mononuclear Cells




**Figure 3. Single Biomarker Cell Mimic Controls Spiked Into PBMCs to Represent Exhausted or Activated States**  
A. Gating scheme showing how cell mimics are represented based on optical parameters on a flow cytometer and gated out from live vs. dead cells.  
B. Representation of Unactivated PBMCs, Activated PBMCs (activated with Dynabeads for 1 week), Unactivated PBMCs were spiked with single biomarker cell mimics expressing CD69 or CD25 to represent activated cells.  
C. Representation of Unactivated PBMCs, Activated PBMCs (activated with Dynabeads for 1 week), Unactivated PBMCs were spiked with single biomarker cell mimics expressing PD-1 or TIM-3 to represent exhausted cells.

## INTRODUCTION

Slingshot leverages the principles of biochemistry, high-precision manufacturing, and polymer chemistry to engineer cell mimics that match the features of biological cells, including optical, fluorescence and biochemical properties for Flow Cytometry. By addressing the limitations of biological controls, cell mimics offer a scalable, reproducible solution to accelerate and improve the development of next-generation immunotherapies.




### Benefits of Cell Mimics



**Figure 1.** Schematic of Slingshot technology capability to match biological cells and benefits for standardized assay development

## Seamless Integration of Slingshot Bioscience's Cell Mimic Controls into your workflow




**Figure 2.** Schematic of Slingshot Bioscience's Cell mimics into Cell therapy workflow

### Advantages:

- **Clone Compatibility:** Fully compatible with commonly used antibody clones across immune profiling panels.
- **Workflow Integration:** Seamlessly integrate into standard staining protocols and multicolor flow cytometry panels.
- **Gating Accuracy:** Match PBMCs in size and internal complexity, supporting existing gating strategies without adjustment.
- **Enhanced Assay Performance:** Improve assay calibration, gating validation, and instrument performance monitoring.

## Development Of Multi-Biomarker Lymphocyte Subset Cell Mimics Demonstrating T, B and NK Cell Populations



**Figure 4.** Schematic of Slingshot Biosciences' capability to match biological cells and representative flow plots for flow cytometry data for TruCytes Lymphocyte Subsets cell mimics compared to cryopreserved PBMCs. All samples underwent identical staining protocols using labeled monoclonal antibodies to CD45, CD3, CD19, CD4, CD8, CD14, CD16, and CD56, were collected on a Cytex Aurora spectral flow cytometer and analyzed via FlowJo.

## Pre Approved Biomarkers for Custom Solutions\*

<b>Activation Markers</b>	CD25, CD69, CD71 (Transferrin Receptor), CD134 (OX40), CD137, CD278 (ICOS), CD80, CD86
<b>Exhaustion Markers</b>	CD152 (CTLA-4), CD279 (PD-1), CD274 (PD-L1), CD223 (LAG-3), TIM-3, TIGIT
<b>Chimeric Antigen Receptors</b>	G4S, Whitlow, Custom CARs
<b>T Cell Markers</b>	CD2, CD3, CD4, CD5, CD7, CD8, CD25, CD27, CD28, CD45, CD154 (CD40L)
<b>T Cell Memory Subset Phenotype markers</b>	CD62L, CCR7, CD45RA, CD45RO, CD27, CD28, CD95
<b>Tumor/Proliferation Markers</b>	Ki-67 (MKI67), CD71 (Transferrin Receptor), CD99, CD227 (MUC1), CD340 (HER2/neu), LRRC15
<b>NK Cell Markers</b>	CD16, CD56, CD335 (NKG46)
<b>Transcription Factors</b>	FOXP3, ROR1
<b>Cytokines &amp; Cytokine Receptors</b>	IL-2, IL-4, IL-21, TNFα, TGFβ2, CD120b (TNF-R1), CD122 (IL-2Rβ), CD124 (IL-4R), CD126 (IL-6R), CD127 (IL-7Ra), CD213a1 (IL-13Ra1), CD213a2 (IL-13Ra2), CD215 (IL-15Ra), CD217 (IL-17R), CD218a (IL-18Ra), CD218b (IL-18Rb)
<b>B Cell Markers</b>	CD19, CD20, CD22, CD23, CD24, CD38, CD79b, CD86, CD267 (TACI), CD268 (BAFF-R), CD269 (BCMA)
<b>Myeloid Markers</b>	CD14, CD16, CD33, CD64, CD66, CD68, CD163, CD172a (SIRPα)
<b>Stem/Progenitor Markers</b>	CD34, CD117 (c-Kit), CD123 (IL-3Ra), CD135 (FLT3)
<b>Adhesion Molecules</b>	CD54 (ICAM-1), CD58 (LFA-3), CD102 (ICAM-2), CD166 (ALCAM)
<b>Other Functional Markers</b>	EGFR, CD36, CD39, CD40, CD47, CD52, CD59, CD93, CD109, CD112, CD147, CD150, CD155, CD162, CD203c, CD221, CD229, CD252 (OX40L), CD254 (RANKL), CD257 (APRIL), CD304 (Neuropilin-1), CD314 (NKG2D), CD319 (SLAMF7), CD328 (Siglec-7), CD329 (Siglec-9)

\*Other markers can be evaluated

**Table 1.** Slingshot evaluated biomarkers for custom projects. These markers can be added to Slingshot cell mimics to develop customizable controls to meet unique assay needs.

## CONCLUSIONS

Slingshot Biosciences has developed precision-engineered cell mimics serve as scalable, reproducible quality controls for flow cytometry, offering several key advantages:

- Bright, specific expression of defined biomarkers with phenotypic and optical profiles that closely resemble activated or exhausted T cells, enabling reliable performance when spiked into PBMC samples.
- Improved assay consistency and reproducibility across a wide range of immune profiling workflows.